

Plasma β -Amyloid 1-42/1-40 Ratio Provides Insight into the Presence of Alzheimer's Disease

Background

- Alzheimer's disease (AD) is a chronic neurodegenerative disease characterized by extracellular brain aggregates of β -amyloid peptides.¹
- These peptides originate from the amyloid precursor protein (APP). Cleavage of APP primarily yields 2 peptides: β -amyloid 1-42 and β -amyloid 1-40.¹
- Measurement of β -amyloid peptide 1-42 in cerebral spinal fluid (CSF), as well as tau and hyperphosphorylated tau, has become a standard diagnostic tool for AD. A blood-based method would avoid the invasiveness of collecting a CSF specimen.
- In addition, some studies suggest that the ratio of 1-42/1-40 in plasma can help identify patients with AD.²⁻⁴
- Objective:** The investigators of this study developed a method to simultaneously measure β -amyloid 1-42 and 1-40 in plasma specimens. They also evaluated the correlation of the 1-42/1-40 ratio with mild cognitive impairment and AD.

Methods

- A highly sensitive multiplex assay for the simultaneous measurement of β -amyloid 1-42 and 1-40 was developed and validated at Quest Diagnostics.
- Clinical specimens were selected from 100 patients who were categorized based on the Clinical Dementia Rating test and the Montreal Cognitive Assessment. The categories were as follows:
 - Normal (n=35)
 - Early mild cognitive impairment (MCI; n=25)
 - Late MCI (n=15)
 - AD (n=25)
- β -amyloid 1-42 and 1-40 were measured using the assay, and 1-42/1-40 ratios were calculated.
- Categories were compared using the Kruskal-Wallis H test.

Results

- β -amyloid 1-42/1-40 ratios were generally higher in specimens from normal patients than those from AD-affected patients ($P=0.008$).
 - Pairing AD with late MCI data and normal with early MCI data slightly improved the distribution ($P=0.001$).
- For identifying patient with AD vs normal patients, the clinical specificity was 71% and sensitivity was 76%.

Conclusions

- The investigators developed an assay that simultaneously measures β -amyloid peptides 1-42 and 1-40 in plasma specimens.
- The assay allows assessment of the 1-42/1-40 ratio, which may help identify patients with AD.

Poster presentation at the Alzheimer's Association International Conference

Authors

Keith R Morneau,¹ Brian G Sansoucy,¹ Robert J Lagier,² Charles M Rowland,² Torey Neusch,¹ Altin Qeleshi,¹ Edward I Ginns,¹ Nigel J Clarke³

Affiliations

¹ Quest Diagnostics, Marlborough, MA
² Quest Diagnostics, Alameda, CA
³ Quest Diagnostics, San Juan Capistrano, CA

Alzheimer's Association International Conference, Los Angeles, CA, July 14-18, 2019

References

- Selkoe DJ, Hardy J. *EMBO Mol Med.* 2016;8:595-608.
- de Rojas I, Romero J, Rodríguez-Gomez O, et al. *Alzheimers Res Ther.* 2018;10:119.
- Fandos N, Pérez-Grijalba V, Pesini P, et al. *Alzheimers Dement (Amst).* 2017;8:179-187.
- Shahpasand-Kroner H, Klafki WH, Bauer C, et al. *Alzheimers Res Ther.* 2018;10:121.